

TP 2x Master Mix

(Cat. No T601, T602, T603, T603xl)

rev. 04/2025

Description

TP 2x Master Mix is dedicated for universal analysis of DNA samples using PCR. It is based on recent finding that addition of Trehalose or 1,2-Propanediol (abbreviated **TP**) into reaction mixture is capable of substantially increasing efficiency of PCR and enhance amplification of samples which are otherwise difficult to amplify, including DNA from whole blood, GC rich amplicons and samples containing PCR inhibitors (*Horáková a spol., BMC Biotechnology, 11:41, 2011*).

Rapid preparation of the samples

- All components of the TP 2x Master Mix are 2x concentrated (optimized reaction buffer containing trehalose and 1,2-propanediol, nucleotides and Taq DNA polymerase) which facilitates rapid preparation of the samples. The samples are prepared by mixing an aliquot of the Master Mix with oligonucleotide primers, template DNA and H₂O (included).
- TP 2x Master Mix is especially useful for routine analyses of large numbers of DNA samples. To 0.5 ml of the Master Mix in original tube, primers (e.g. 40 µl forward and 40 µl reverse) and PCR H₂O (e.g. 380 µl) are added and mixed; the "armed" Mix can be stored at -20 ± 5°C. Immediately before use, the Mix is thawed and each e.g. 24 µl aliquot is mixed with e.g. 1 µl of the tested DNA template.

Technical data

Components and packaging

- 1 tube with 0.5 ml TP 2x Master Mix (for 40 reactions, 25 µl each)
- 1 tube with 1.5 ml PCR H₂O.

Storage

- For very short terms (hours, days) at room temperature
- For short terms (weeks) at 4 - 8°C.
- For long terms at -20 ± 5°C. Material can be repeatedly defrosted

Composition

- TP 2x Master Mix contains: 150 mM Tris-HCl, pH 8.8 (25°C), 40 mM (NH₄)₂SO₄, 0.4 M trehalose, 2 M 1,2-propanediol, 0.02% Tween 20, 5 mM MgCl₂, 400 µM dATP, 400 µM dCTP, 400 µM dGTP, 400 µM dTTP, Taq DNA polymerase (50 U/ml), stabilizers and additives.

Purity and quality control

- Purity of Taq DNA polymerase is verified by SDS PAGE, only one band of 94 kDa is observed in Coomassie blue stained gel. Material is free of nucleases.
- Each batch of TP 2x Master Mix is tested for amplification of a single copy gene with high content of GC in genomic DNA.

Cat. No.	Product name	Quantity
T601	TP 2x Master Mix (1x)	40 reactions
T602	TP 2x Master Mix (5x)	200 reactions
T603	TP 2x Master Mix (25x)	1000 reactions
T603xl	TP 2x Master Mix (100x)	4x 1000 reactions



Protocol

Suggested basic protocol for PCR amplification using TP 2x Master Mix

1. In a thin-walled PCR tube the following components are mixed

Component	PCR in 25 μ l*	Final concentration
TP 2x Master Mix**	12.5 μ l	75 mM Tris-HCl, pH 8.8, 20 mM (NH ₄) ₂ SO ₄ , 0.2 M trehalose, 1 M 1,2-propanediol, 0.01% Tween 20, 2.5 mM MgCl ₂ , 200 μ M each of dNTPs, 25 U/ml Taq DNA polymerase, stabilizers and additives
5' primer (50 μ M)	1 μ l	2 μ M
3' primer (50 μ M)	1 μ l	2 μ M
Template DNA (1 ng/ μ l - 1 μ g/ μ l); or nonseparated 2x diluted blood	1 μ l	0.04 ng – 0.04 μ g DNA/ μ l
PCR H ₂ O (Cat. No. P042)	9.5 μ l	

* Different volumes can be used, but the Master Mix should be finally diluted twice.

** Before using frozen Master Mix it is important that all components are completely soluble. Solubilization is accelerated by warming Master Mix at 37°C.

2. Mix gently and briefly centrifuge.

3. Perform PCR under conditions optimized for the primers used. Common cycling parameters are:

I. Initial denaturation, 94°C, 10 min

II. Cycling and amplification of the template

Denaturation 94°C, 10 sec

Primers annealing 55 - 65°C (depending on the primers), 10 sec

Extension 72°C, 10-30 sec (~20 sec for 500 bps)

During this step fluorescence of SYBR I is measured

Repeat 30 – 40x

III. Amplified DNA can be loaded into agarose gel after adding a loading buffer (Cat. No P048, P062, P064 or P066).